

# GAMMEX® Non-Latex PI Micro

S U R G I C A L   G L O V E S

## Polyisoprene, skin-friendly surgical glove with exceptional dexterity and tactile sensitivity for surgical finesse

### KEY FEATURES AND BENEFITS:

- 15% thinner on average than Ansell's standard polyisoprene surgical gloves for unsurpassed dexterity and tactile sensitivity
- Softer and 'stretchier' PI glove featuring superior comfort compared to Ansell's latex gloves
- Smooth with micro-textured finish for exceptional instrument handling and grip

### Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

### FEATURED TECHNOLOGIES



## GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Non-Latex PI Micro

SURGICAL GLOVES

## PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	White
Thickness	Micro
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D7160, EN 16523-1, EN 374:2003, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340057055	340057060	340057065	340057070	340057075	340057080	340057085	340057090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu

**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com