

MICRO-TOUCH[®] SensiClean[®] II

Essential Series

Latex, powder-free, examination glove with heightened tactile sensitivity and dexterity

This non-sterile, glove combines excellent tactile sensitivity with high levels of comfort while delivering a positive, non-aggressive grip. The natural rubber latex formulation ensures a comfortable, well-fitting glove while Ansell's extensive leaching phases ensure a very low residual protein content, therefore suitable for sensitive, frail skin.

KEY FEATURES AND BENEFITS:

- Double chlorinated inner surface is ideal for effortless double gloving while remaining very gentle on fragile skin
- Mildly textured surface for secure wet and dry grip
- Low allergenicity latex reduces risk of contact dermatitis due to absence of thiurams

Caution: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.

Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Brief, low-risk examination procedures
- Damp donning



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] SensiClean[®] II

Essential Series

PRODUCT INFORMATION

| | | | |
|---|----------------------|------------------------------------|--|
| Material | Natural Rubber Latex | Not Made From Natural Rubber Latex | No |
| Color | White | Cuff Length | Standard |
| Powder Content | Powder-Free | External Glove Surface | Fully Textured |
| Freedom from Holes (Inspection level I) | 1.5 AQL | Palm Thickness (mm/mil) | 0.14 / 5.5 |
| Finger Thickness (mm/mil) | 0.15 / 5.9 | Available Sizes | XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10) |
| Tested For Use With Chemotherapy Drugs | No | Sterile | No |
| Antistatic | Not Tested | Glove Length (mm/inches) | 240 / 9.5 |
| Product Segmentation | Essential Series | | |

ORDERING INFORMATION

| Size | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | 4566 | 4560 | 4562 | 4564 | 4568 |

PACKAGING AND STORAGE

| | |
|----------------------|--|
| Packaging | 100 gloves per dispenser 10 dispensers per case 1000 gloves per case |
| Shelf Life | 3 Years |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

STANDARDS AND CERTIFICATIONS

EN 455 Part 2, EN 455 Part 3, EN 455 Part 1, EN 420:2003 + A1:2009, ISO 11193, ISO 9001, ISO 13485, ISO 10993-10

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com