

Transparent Silicone Mask



Features

- Fit snugly around the face to minimize the dead space.
- High transparency for patient's visibility.
- Reusable 1- piece mask in Silicone offers soft and optimal seal.
- Autoclavable up to 134°C.

Code	Description	Package
SM-86583	Silicone Mask, Size 000	150/case
SM-86573	Silicone Mask, Size 00	150/case
PN-0000	Silicone Mask, Size 0	300/case
PN-0001	Silicone Mask, Size 1	150/case
PN-0002	Silicone Mask, Size 2	125/case



Features

- Liquid injection assures steady color without distortion after sterilization.
- Reusable 1- piece mask in Silicone offers soft and optimal seal.
- High transparency for patient's visibility.
- Autoclavable 134°C.

Code	Description	Package
PN-0006A	Silicone Mask, Size 3, Child Large	125/case
PN-0007A	Silicone Mask, Size 4, Adult	125/case
PN-0008A	Silicone Mask, Size 5, Adult Large	50/case



Silicone Mask

Features

- The masks are designed to fit closely around the face to minimize dead space.
- 100% Medical Silicone, Polycarbonate, Latex free.
- Autoclavable up to 121°C.

Code	Description	Package
PN-0003	Silicone Mask, Size 3, Child Large	150/case
PN-0004	Silicone Mask, Size 4, Adult	100/case
PN-0005	Silicone Mask, Size 5, Adult Large	100/case

