



**Australian Government**  
**Department of Health and Aged Care**  
Therapeutic Goods Administration

## Public Summary

**Summary for ARTG Entry:** 157331 Convatec Australia Pty Ltd - Duoderm Extra Thin - Dressing, hydrogel

**ARTG entry for** Medical Device Included Class III  
**Sponsor** ConvaTec Australia Pty Ltd  
**Postal Address** PO Box 63, MULGRAVE, VIC, 3170  
Australia  
**ARTG Start Date** 28/11/2008  
**Product Category** Medical Device Class III  
**Status** Active  
**Approval Area** Medical Devices

### Conditions

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.  
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

### Manufacturers

Name	Address
Convatec Limited	Unit 20 1st Avenue Deeside Industrial Park Deeside, Flintshire, CH5 2NU United Kingdom

### Products

#### 1 . Duoderm Extra Thin - Dressing, hydrogel

<b>Product Type</b>	Single Device Product	<b>Effective Date</b>	28/11/2008
<b>GMDN</b>	34082 Dressing, hydrogel		
<b>Functional Description</b>	Duoderm Extra Thin dressings interact with wound moisture producing a soft mass that enables removal of the dressing with little or no damage to newly formed tissues. They help isolate the wound against bacteria and other external contamination		
<b>Intended Purpose</b>	Duoderm Extra Thin is a highly flexible, control gel formula dressing designed for use on dry to lightly exudating wounds.		
<b>Variant information</b>	Size (cm) .4-35cm		

### Specific Conditions

No Specific Conditions included on Record

© Commonwealth of Australia. This work is copyright. You are not permitted to re-transmit, distribute or commercialise the material without obtaining prior written approval from the Commonwealth. Further details can be found at <http://www.tga.gov.au/about/website-copyright.htm>.

Public Summary