

HipSaver – the original soft hip protector with proven protection

HipSaver Selection Guide

Wearer's Particulars	Nursing Home Undergarment	SlimFit Undergarment	QuickChange Undergarment	Open-Bottom Undergarment	Interim Overpant	Hip Protecting Track Pant	Hip Protecting Shorts
Socially active and is cognitively well		✓				✓	✓
Prefers a discreet garment		✓					
Continent		✓				✓	✓
Wears an incontinence pad	✓	✓	✓		✓	✓	✓
Wears a bulky incontinence brief/diaper	✓		✓	✓	✓	✓	✓
Experiences urge incontinence				✓			
Staff changes incontinence product in bed			✓				
Tends to fall forward - Knee padding optional						✓	
Falls backwards - Tailbone padding optional	✓	✓	✓		✓	✓	✓
Wearer specific due to hygiene regulations	✓	✓	✓	✓			
Multi-wearer interchangeable after laundering					✓	✓	✓
Immediate protection for new admissions					✓	✓	✓

HipSaver Model Descriptions

Pant Only: 1 pant with open pouches for holding removable interchangeable hip protecting pads.
(hip protecting pads not included)

All Soft Hip Protecting Pads: 1 set of hip protecting pads with AirPad technology. (i.e. 2 pads)

Starter Kit: 3 pants with open pouches and 1 set of removable interchangeable hip protecting pads.

Veterans Kit 2: 6 pants with open pouches and 2 sets of removable interchangeable hip protecting pads.

High Compliance: 1 pant with hip protecting pads permanently sewn in.

High Compliance with Tailbone Protection: 1 pant with hip and tailbone protecting pads permanently sewn in.

Veterans High Compliance Pack: 6 pants with hip protecting pads permanently sewn in.

All HipSaver products including soft hip protecting pads are fully machine washable and tumble dryable in both home and institutional laundries up to 95°C.

HipSaver Measurement Chart

Measuring Instructions	Hip Measurement		Size
	cm	inches	
Nursing Home, SlimFit, and QuickChange styles are available in 'Male' with a flyfront or 'Female' without. Please specify gender when ordering these styles. All other garments are unisex. When assessing wearer's size, measure all the way around the widest part of the hips, usually over the greater trochanters. Measure over any incontinence product if applicable. Please measure carefully as hygiene regulations prevent us from accepting returned garments once they have been worn. Manufacturing anomalies excepted.	71 to 81 cm	28 to 31 inches	Extra Small
	82 to 92 cm	32 to 35 inches	Small
	93 to 102 cm	36 to 39 inches	Medium
	102 to 114 cm	40 to 44 inches	Large
	115 to 127 cm	45 to 50 inches	Extra Large
	128 to 146 cm	51 to 57 inches	Extra Extra Large