

# VBag®

# Sickeze®

## VBAG® EMESIS BAG & SICKZEZE® RETAIL SHELF PACK



Sentry Medical Emesis Bags have a firm polypropylene twist and lock rim to provide easy control, and help prevent spillage after use. They isolate splash back and have a graduated indicative measurement system to 1500ml, allowing measurement of contents in the bag. On completion of use, the bag is easily disposed of with no mess. Simply twist the bag and lock into the keyhole slot on the rim.

### PRODUCT FEATURES

- o Rims are colour coded red for infection control
- o Convenient dispenser also available
- o Easy twist and lock to seal
- o Assists in infection control

### SUITABLE FOR

- o Infection controlled area
- o First aid
- o Hospital and medical



REORDER	PRODUCT	PCS/UNIT	UNITS/CTN
VB002	VBag® Emesis Bag	50	10
VB011	Sickeze® Emesis Bag Retail Shelf Pack (3pcs/pk - 6pks)	10	1
VBD001	Acrylic Emesis Bag Dispenser	1	12

All product information is correct at date of publication

