



STAINLESS STEEL HOSPITAL UTENSIL & MEDICAL WARE

Diverse surgical stainless steel products for patient care

The Liberty® branded stainless steel ware is made from a high quality genuine stainless steel 18-8 (304) specifically manufactured for surgical products with seamless construction that prevents bacteria growth and allows products to be cleaned with ease.

Manufacturing premium grade stainless steel hospital utensils and medical equipment that meets the most exacting requirements of medical and healthcare professionals worldwide. The medical grade stainless steel is a highly corrosion and rust-resistant alloy. The metal is extremely strong and durable and, more importantly, has a characteristic ability to form protective layers.



The Liberty® range of stainless steel hospital holloware includes steel trays for different purposes, buckets, lids, kidney dishes, bowls, cups, jugs, jars, instrument holders and bed pans, etc. are preferred in steel. These products are manufactured using seamless technology which eliminates areas that can harbour bacteria growth.



Code	Description	Size
8002	Bucket Stainless Steel	300 x 265mm
8003	Bucket Stainless Steel	265 x 310mm
8003CP	Lid for 8003 (not pictured)	

Liberty®



STAINLESS STEEL HOSPITAL UTENSIL & MEDICAL WARE



Code	Description	Size	Code	Description	Size
K037	Instrument Tray only	340 x 240 x 50mm	M013	Instrument Tray with Lid	300 x 200 x 50mm
K050	Instrument Tray Shallow	480 x 330 x 19mm	M018	Kidney Dish	172 x 93 x 34mm
M002	Catheter Tray with Lid	200 x 80 x 38mm	M018L	Kidney Dish Lid only	177mm
M003	Kidney Dish	200 x 95 x 38mm	M019	Catheter Tray with Lid	430 x 100 x 50mm
M003L	Kidney Dish Lid only	210mm	M028	Instrument Tray with Lid	300 x 200 x 65mm
M005	Instrument Tray with Lid	200 x 130 x 50mm	M034	Kidney Dish	300 x 145 x 57mm
M006	Kidney Dish	250 x 110 x 45mm			
M006L	Kidney Dish Lid only	255mm			



STAINLESS STEEL HOSPITAL UTENSIL & MEDICAL WARE



Code	Description	Size
K082	Lotion Bowl	140 x 60mm
K083	Lotion Bowl	160 x 70mm
K084	Lotion Bowl	180 x 80mm
K125	Foot Basin (not pictured)	350 x 115mm
K128	Wash Bowl 3.7 Litre	320 x 85mm
M042	Iodine Bowl 55cc	65 x 33mm
M043	Iodine Bowl 110cc	77 x 43mm

Code	Description	Size
M044	Iodine Bowl 170cc	83 x 50mm
M045	Iodine Bowl 220cc	92 x 52mm
M046	Iodine Bowl 280cc	100 x 52mm
M048	Iodine Bowl 390cc	120 x 60mm



STAINLESS STEEL HOSPITAL UTENSIL & MEDICAL WARE



Code	Description	Size
M001	Medicine Cup	2oz
M009	Ointment Jar with Lid	75 x 57mm
M023	Bed Pan	Adult
M024	Douche Can	2 Litre
M025	Graduated Jug	2000cc
M026	Dressing Jar with Lid	75 x 57mm
M029	Dressing Jar with Lid	100 x 125mm
M030	Graduated Jug	1000cc

Code	Description	Size
M031	Forceps Jar	55 x 114mm
M032	Forceps Jar	55 x 190mm
M033	Graduated Jug	500cc
M038	Thermometer Jar	25 x 100mm