

## Fingertip Pulse Oximeter

Model: SB100



- %SpO2 & Pulse Rate - Big LED screen display for easy reading
- Shielded Design blocks Ambient Light
- Biocompatibility & Anti-Allergic Design
- One- piece face cover design & One-touch ON/OFF key
- Pulse Strength Indicator
- Lanyard-attached - Economical & portable health management tool



%SpO2 & Pulse Rate



Large Font Size



Shielded Design blocks Ambient Light



Biocompatibility & Anti-Allergic Design



Pulse Strength Indication



String-Attached

### Specification

**Measurement Range** | SpO2: 35%~99%; PR: 30~250 bpm

**Precision** | SpO2: 70%~99%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3

**Power Supply** | 2 x "AAA" Alkaline

### Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

Model	Qty per carton	Carton volume
SB100	40 pcs	0.028 cbm/ctn

- Storage temperature: -25°C - 70°C (-13°F - 158°F)
- Relative humidity: 15% - 90% (no condensing)
- Atmospheric Pressure: 700 to 1060 hPa