



Water Resistant



Contours to Body



Transparent Film

A transparent adhesive polyurethane film with moisture vapour permeability which adheres to the skin and forms a flexible and comfortable layer. AsGuard Clear conforms to the skin and allows easy visibility to wound healing.

## CUT TO LENGTH

Measured grid backing allows cut to correct size as required

## TRANSPARENT FILM

Allows visualisation of covered tissue

## SHOWER AND INFECTION RESISTANT

Impermeable to help protect site against infection and allows showering with the dressing in-situ



## Suitable for:

Wounds with little or no exudate

Protection and securement of IV lines

Dressing reinforcement

Reorder	Size	Pcs/Unit	Units/Ctn
ISB010	5cm x 10m	1	24
ISB011	10cm x 10m	1	12

