



Strong and Durable



Various Sizes



Colour Change Indicator

Produced with 60gsm medical grade paper and a PET/PP blend for the film, Bolsaplast Sterilisation Reels offer excellent peel ability and is resistant to film tearing. The reels are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation process.

## COLOUR CHANGE INDICATOR

Chemical indicators for Steam / EO/ Formaldehyde

## NON-TOXIC INK

No damage to instruments

## MEDICAL GRADE CONSTRUCTION

Offers excellent peel ability and is resistant to film tearing



## Suitable for:

Sterilising devices and equipment

Reorder	Size	Pcs/Unit	Units/Ctn
MRP1040	5cm x 200m	1	6
MRP1041	7.5cm x 200m	1	4
MRP1042	10cm x 200m	1	4
MRP1043 <<	12.5cm x 200m	1	4
MRP1044	15cm x 200m	1	4
MRP1045	20cm x 200m	1	2
MRP1046	25cm x 200m	1	2
MRP1047	30cm x 200m	1	2
MRP1048	40cm x 200m	1	1

<< Indent Item