



Reusable



Comfortable



Non-Slip

Lymphlex Air is a fine pored, soft-to-touch reusable foam bandage for use as padding under compression bandages and can be washed up to 50 times, making it a cost effective component in oedema treatment.

Lymphlex Air is an individually replenishable component of the Lymphlex Reduce Compression Therapy Kit.

SOFT AND CONFORMING

Provides comfort and protection

NON-SLIP FABRIC

Durable and tear resistant which provides even distribution of pressure

COST EFFECTIVE

Washable and reusable up to 50 times



Suitable for:

Fragile skin

Padding for compression bandage therapy

Aids in helping to reduce oedema in Lymphodema treatment*

Increasing tissue pressure to avert oedema fluid occurrences*

*When used with Lymphlex compression bandages

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
LOR021	8cm x 2.5m	Non-Sterile	1	24
LOR022	10cm x 2.5m	Non-Sterile	1	24

