



Comfortable
for Patient



Fast
Absorbency



Gentle
on Skin

Mediam 210 Low Adherent Dressing consists of an absorbent cotton inner pad encased between two layers of high density polyethylene film.

LOW-ADHERENT

Low adherence to the wound, therefore minimal trauma and pain upon removal

HIGHLY ABSORBENT

Highly absorbent inner pad, wicks excess exudate away from wound bed

CUT TO SIZE

Maintains integrity when cut to size

Suitable for:

Light to moderate exuding wounds

Cuts and grazes

Minor burns

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
NAD001	5cm x 5cm	Sterile	50	10
NAD002	7.5cm x 7.5cm	Sterile	50	10
NAD003	10cm x 10cm	Sterile	50	8
NAD004	10cm x 20cm	Sterile	50	8

