

# NEW 1L Wall Dispensers

Product specifications and dimensions  
(For use with MICROSHIELD® 1L products\*)



## 1L Wall Dispenser - Manual

<b>Product Code</b>	70002053
<b>GTIN</b>	9350279000913
<b>Dimensions</b>	(w) 108mm x (d) 100mm x (h) 310mm
<b>Weight</b>	450 grams
<b>Sales unit</b>	Each
<b>Cleaning and maintenance:</b> Clean exterior surface of dispenser with warm water and mild detergent. DO NOT USE ABRASIVE CLEANERS. Dry with a soft cloth.	



## 1L Wall Dispenser - Automatic/touch-free

<b>Product Code</b>	70002052
<b>GTIN</b>	9350279000906
<b>Dimensions</b>	(w) 118mm x (d) 100mm x (h) 336mm
<b>Weight</b>	700 grams
<b>Batteries</b>	4 x D ( <i>recommend industrial or heavy duty</i> )
<b>Sales unit</b>	Each



## Drip Tray for 1L Wall Dispenser

<b>Product Code</b>	60000379
<b>GTIN</b>	9350279000944
<b>Dimensions</b>	(w) 110mm x (d) 100mm x (h) 135mm
<b>Weight</b>	45 grams
<b>Sales unit</b>	Each

\* Suitable for use with MICROSHIELD ANGEL CLEAR® Antimicrobial Hand Gel 1L and MICROSHIELD® Antiseptic Hand Rub 1L pack