



Highly Absorbent



Low Linting



Up To Moderate Exudate

Non-Woven Combine Dressings consist of a highly absorbent inner, enclosed in a soft and comfortable non-woven covering fabric with soft end seals. Ideal for use as a primary or secondary dressing and are non-abrasive on the skin.

HIGHLY ABSORBENT

Suitable for up to moderate exudate and aids in less dressing changes

ABSORBENT INNER PAD

Inner pad wicks excess exudate away from wound bed



NON-ABRASIVE OUTER

Suitable for fragile and sensitive skin

Suitable for:

Primary or secondary dressing

Moderate exudating wounds

First aid

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
NWC001	10cm x 10cm	Sterile	50	16
NWC002	10cm x 20cm	Sterile	50	8
NWC003	20cm x 20cm	Sterile	50	5

