



Low Linting



100% Cotton



X-Ray Detectable

ThreadX Standard and Optimum Gauze Swabs are 100% cotton, with various weave count options for preferred absorbency. A black x-ray detectable thread heat bonded to the gauze.

### FRENCH FOLD

Helps reduce edge fraying

### X-RAY DETECTABLE

X-ray detectable thread for patient safety

### LOW LINTING

Aids in infection control and minimises the risk of contamination



## Suitable for:

All open surgical procedures

Operating room

Reorder	Ply	Size	Weave Count	Pcs	Sterility	Pcs/Unit	Units/Ctn
<b>XRD020</b>	8ply	10cm x 10cm	Optimum	5pcs	Non-Sterile	600	1
<b>XRD3252</b>	8ply	10cm x 10cm	Standard	5pcs	Sterile	250	1
<b>XRD3253</b>	8ply	10cm x 10cm	Optimum	5pcs	Sterile	250	1