

Endoscopy Water Quality Requirements for Compliance AS/NZS4187:2014

In order to assist endoscopy departments work towards compliance, Whiteley Medical is offering water testing kits and endotoxin compliance testing.

Monthly endoscopy water testing for compliance with AS/NZS4187:2014
Monthly testing is required for water hardness, conductivity and chloride in endoscopy units.

Water Hardness Test Kit

Product Code: 303004

\$130 + GST



HACH Chloride Test Strips (40 strips)

Product Code: 303022

\$191 + GST



*Suitable for testing Chloride levels between 30-600ppm.

Conductivity Meter

Product Code: 303025

\$296 + GST



Table 7.2: AS/NZS4187:2014 WATER QUALITY USED FOR PROCESSING RMDs		
Substance	Water sample	Maximum concentration levels
Water Hardness	Supply Water	150mg/L
Chloride	Supply Water	120mg/L
Conductivity	Final Rinse Water	<30µs/cm



Please scan the QR code for more information on water testing

Annual endoscopy endotoxin testing for compliance with AS/NZS4187:2014			
Final Rinse Water Endotoxin Testing*	Endotoxin (<30 EU/mL)	Annual Endotoxin Testing Fee Product Code: 303024	\$200 + GST

*Microbiological Compliance testing is also required monthly and should be conducted with your normal microbiology test provider.

For more information contact whiteley@whiteley.com.au or 1800 833 566

Insist on Whiteley Medical for your Water Testing information

For more information about Whiteley Medical products visit www.whiteley.com.au or call the Product Support Hotline on **1800 833 566**

© Registered Trademark Whiteley Corporation Pty. Ltd. © 2019

