



AAMI  
Level 2



Various  
Colours



Impervious

Isolation Gowns comply with ISO16603, ISO16604 and ISO6529 standards, as well as meet AAMI Level 2 liquid barrier performance levels.

**IMPERVIOUS**

An impervious layer prevents fluid to pass through

**SOF-CUFF**


Added comfort and ease when donning gloves

**VELCRO NECK AND TIE WAIST**


Velcro neckline and back ties to easily don and doff



Blue



Yellow



Purple  
(Oncology)

**Suitable for:**

Hospital and medical

Laboratory

Reorder	Size	Feature	Colour	Sterility	Pcs/Unit	Units/Ctn
GN001P	Regular	AAMI Level 2	Blue	Non-Sterile	50	1
GN003P	Regular	AAMI Level 2	Yellow	Non-Sterile	50	1
GN0041P	Regular	AAMI Level 2	Purple	Non-Sterile	50	1