

## Underpad - 4 ply



**SOFT**



**ABORBENT**



**DRY**

<b>ITEM DESCRIPTION:</b>	Underpad 4ply Blue
<b>SIZE:</b>	40cm x 60cm
<b>SELL UOM:</b>	1 ctn 300 (12pkts of 25)
<b>PRODUCT CODE/S:</b>	CCUP4060
<b>ARTG#:</b>	383870
<b>MSDS:</b>	Available on request
<b>MATERIAL:</b>	1st ply: Non-woven 12gsm; 2nd ply: Paper tissue 14gsm; 3rd ply: Import Fluff pulp; 4th ply: Paper Tissue 14gsm; Cover: 22gsm blue extruded PE film.
<b>WATER ABSORPTION:</b>	110~160ml.
<b>STERILITY:</b>	Non-Sterile



The underpads provides soft, non-woven fabric with a highly absorbent fill material to maximise ability to absorb unwanted fluids. The underpads contains an absorbent inner to lock away wetness. It is sealed on four sides, with a waterproof backing to prevent strikethrough.

- Absorbent – Absorbs and holds in liquid.
- Soft top cover fabric to increase comfort while maintaining its strength
- Soft waterproof backing is silent and minimise movement when in use.
- Waterproof backing prevents leaking

#### SUITABLE FOR:

- Aged Care
- Baby changing mat
- Chair and bed protection
- Incontinence Protection

#### MARKETS:

- General Practice
- Aged Care
- Hospital

