

Disposable Vinyl Exam Gloves



CLEAR

BLUE

FACTS AND FEATURES:

- Powder free
- Latex free formulation, no natural rubber protein.
- Chemical accelerator free.
- Non Powdered Interior Polyurethane coated for easing donning.
- Smooth finish for tactile sensitivity

GLOVE ATTRIBUTES

Design	Ambidextrous, non sterile, beaded cuff
Donning agent	Non Powdered Polyurethane coating on interior surface.
Colour	Clear, Blue
Selling UOM	Box/100
Usage	Single use only
ARTG#	376816

CODES

SIZE	CLEAR VINYL EXAM GLOVES	BLUE VINYL EXAM GLOVES
Small	EBVGCS100	EBVGBS100
Medium	EBVGCM100	EBVGMB100
Large	EBVGCL100	EBVGBL100
X-Large	EBVGCXL100	EBVGBXL100

GLOVE PROFILE

VINYL EXAM GLOVES

Length (mm)	Minimum 230 mm		
Thickness (mm)	Cuff (mm) :Minimum 0.05 mm	Palm (mm): Minimum 0.08 mm	Finger (mm): Minimum 0.08 mm
Tensile Strength(MPa)	Before aging: Minimum 11MPa	After aging: Minimum 11MPa	
Elongation(%)	Before aging: Minimum 350%	After aging: Minimum 350%	
Width (+/- 5 mm)	S= 85 mm M = 95 mm L= 105 mm XL = 115 mm		
Weight (+/- 0.3 gm)	Small = 5.0gm Medium = 5.5gm Large = 6.1gm X-Large = 6.9 gm		



Distributed By:



www.eboshealthcare.com.au