

**NULINE**  
REFRIGERATION

## Win a Free Medical Fridge or Freezer\*



We are searching for the oldest working Nuline Refrigeration's medical fridge or freezer as part of our 40<sup>th</sup> Anniversary celebration.

Win a new Nuline Refrigeration's medical fridge most like your existing one.

**Proudly Manufacturing in Australia and a locally owned family business.**

Simply fill in the form on our website

Submissions close June 30<sup>th</sup> 2025, and the winning entrant will be notified via email and phone on July 14<sup>th</sup> 2025.

\* Terms and Conditions apply



# NLM SERIES LABORATORY REFRIGERATOR

**NULINE**  
REFRIGERATION

## FEATURES

- Top Mounted System
- E-Glass and Heated Cabinet Facia to prevent condensation up to 80% R.H.
- Interior Lighting LED
- Automatic Water Disposal
- Digital Display with Built-in Alarms
- ARTG Number 269894



2°C to 8°C



400 L  
700 L  
1000 L  
1614 L  
2170 L



600w x 620d x 1915h  
870w x 700d x 2100h  
1200w x 700d x 2100h  
1880w x 700d x 2100h  
2400w x 700d x 2100h

+61 3 9562 3355  
info@nuline.au  
[www.nuline.au](http://www.nuline.au)



# NLM SERIES LABORATORY REFRIGERATOR



Model	NLM400/1	NLM700/1	NLM1000/2	NLM1614/3	NLM2170/4
Gross Capacity	400 L	700 L	1000 L	1614 L	2170 L
External Dimensions in mm (w/d/h including castors)	600 x 620 x 1915	870 x 700 x 2100	1200 x 700 x 2100	1880 x 700 x 2100	2400 x 700 x 2100
Internal Dimensions in mm (w/d/h)	500 x 470 x 1435	770 x 550 x 1555	1100 x 550 x 1555	1780 x 550 x 1520	2300 x 550 x 1520
Packed Dimensions in mm (w/d/h including pallet)	700 x 720 x 2065	970 x 800 x 2250	1300 x 800 x 2250	1980 x 800 x 2250	2500 x 800 x 2250
Gross Weight - kg (packed weight)	97	107	137	187	207
Net Weight - kg	85	95	125	175	195
Voltage	240V				
Temperature Range	2°C to 8°C				
Ambient Temperature Range	10°C to 32°C				
Refrigerant	R134a				
Cooling System / Defrost	Fan Forced / Automatic				
Exterior Cabinet Finish	White Colourbond				
Inner Liner Material	White Colourbond				
Door	Solid / Glass E-Glass				
Door Hinge	RH		LH / RH	LH / RH / RH	LH / RH / LH / RH
Door Handle	Lockable				
Type of Controller	Digital				
Temperature Display	External				
Data Logger	Optional				
Door Ajar Alarms	Optional				
High and Low Readout	Visual / Audible				
BMS Connectivity					
(BMS Connectivity for Remote Alarms only)	Optional				
(BMS Connectivity for 4-20mA Transmitter)	Optional				
Port Entry	Optional				
Battery Back-up for Controller	Optional				
Power Requirement	10amp GPO				
Castors	Yes				
Adjustable Shelves		5		15	20
Usable Shelf Area in mm(w/d)	480 x 450	750 x 530	1080 x 530	565 x 530	575 x 530
Light	LED				

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

# AUSTRALIAN MADE

The logo for Nuline Refrigeration is contained within a blue oval with a white border. The word "NULINE" is written in a large, white, sans-serif font, and the word "REFRIGERATION" is written in a smaller, white, sans-serif font directly below it.

**NULINE**  
REFRIGERATION

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to  $-40^{\circ}\text{C}$ ), blanket warmer, mortuary refrigerators and refrigerated incubators.

