



Strong and Durable



Various Sizes



Colour Change Indicator

Produced with medical grade paper on one side and a PET/PP blend film on the other, Bolsaplast Sterilisation Gusset Reels offer extremely low fibre shedding and are resistant to film tearing. The reels are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.

## COLOUR CHANGE INDICATOR


Chemical indicators for Steam/EO/Formaldehyde

## NON-TOXIC INK

No damage to instruments

## MEDICAL GRADE CONSTRUCTION

Offers excellent peelability and is resistant to film tearing



## Suitable for:

Sterilising devices and equipment

Reorder	Size	Pcs/Unit	Units/Ctn
MRF1060	10cm x 4cm x 100m	1	4
MRF1061	15cm x 4cm x 100m	1	4
MRF1062	20cm x 6cm x 100m	1	2
MRF1063	25cm x 6cm x 100m	1	2
MRF1064<<	30cm x 8cm x 100m	1	1
MRF1065<<	40cm x 8cm x 100m	1	1

<< Indent Item