

# GAMMEX<sup>®</sup> Latex

S U R G I C A L   G L O V E S

## Latex, powder-free surgical gloves combining optimal comfort, strength and tactile sensitivity

- **Elevated comfort and protection:** GAMMEX<sup>®</sup> Latex powder-free gloves' unique latex formulation is soft yet robust; delivering the best combination of comfort and hand protection
- **Certified for ergonomic comfort and performance:** Certified by United States Ergonomics, this gloves' anatomically designed formers feature finger shaped pads that heighten fingertip sensitivity resulting in minimal hand muscle effort used
- **Extra cleaning processes:** Additional cleaning processes ensure the lowest levels of latex proteins and chemical accelerators to help minimize the risk of developing Type I latex or Type IV chemical allergies
- **Effortless donning:** Ansell's unique donning aid technology helps make dry and damp hand donning effortless
- **Safe chemo handling:** Tested in accordance with ASTM D6978 and may be safely used across a range of chemotherapy drugs\*



**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

### Recommended For

- Damp donning
- Extra comfort during long, sustained procedures
- Tested for use with chemotherapy drugs
- Use as an outer glove when double gloving

### FEATURED TECHNOLOGIES



## GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex

## SURGICAL GLOVES

### PRODUCT INFORMATION

|                                                |                                                                                                                                                            |
|------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>                                | Natural Rubber Latex                                                                                                                                       |
| <b>Powder Content</b>                          | Powder-Free                                                                                                                                                |
| <b>Color</b>                                   | Natural                                                                                                                                                    |
| <b>Thickness</b>                               | Standard                                                                                                                                                   |
| <b>Cuff Style</b>                              | Straight with SUREFIT™ Technology                                                                                                                          |
| <b>External Glove Surface</b>                  | Smooth with Micro-Textured Finish                                                                                                                          |
| <b>Internal Glove Surface</b>                  | Polymer Coated with DERMASHIELD™ Technology                                                                                                                |
| <b>Allergy Prevention</b>                      | None                                                                                                                                                       |
| <b>Grip Level</b>                              | Moderate                                                                                                                                                   |
| <b>Double Gloving Recommendation</b>           | Outer glove                                                                                                                                                |
| <b>Tested For Use With Chemotherapy Drugs</b>  | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)                                                                                         |
| <b>Freedom from Holes (Inspection level I)</b> | 0.65 AQL                                                                                                                                                   |
| <b>Sterile</b>                                 | Yes                                                                                                                                                        |
| <b>Product Standards</b>                       | AS/NZS 4179, ASTM D6978, ASTM F739, EN 16523-1, EN 374-2 and -4, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107 |
| <b>Quality / Environmental Standards</b>       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP                                                                                       |
| <b>Packaging</b>                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                                                                       |
| <b>Shelf Life</b>                              | 3 years                                                                                                                                                    |
| <b>Storage Instructions</b>                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                                                   |
| <b>Vulcanization Chemical Accelerators</b>     | Xanthogen Polysulphides, Sodium di-n-butylthiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC)                                                         |

### ORDERING INFORMATION

| Size                | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         | 9.5       |
|---------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| <b>Product Code</b> | 330048055 | 330048060 | 330048065 | 330048070 | 330048075 | 330048080 | 330048085 | 330048090 | 330048095 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com